

RESTING ECG

SIGNAL QUALITY AS A KEY FOR PRECISE EVALUATION

All the BTL electrocardiographs were designed especially for simple operation, suitable even for the most demanding experts.

In combination with BTL CardioPoint-ECG software, the BTL ECGs allow comparing ECG records, sharing them and exporting them to third-party medical systems.

POWERED BY



PROFESSIONAL EXTENTIONS



Automatic ECG measurements



Spirometry



WiFi module



User-friendly software



Sudden death screening



Data sharing


BTL CARDIPOINT ECG



BTLNET.COM



CHOOSE WHAT FITS YOU



	L-LINE ECG	M-LINE ECG	S-LINE ECG
Display	Colour touch 8.4" / 5.7"	Colour touch 5.7"	Graphical
Printing paper	A4/A5 Z-fold, roll, fax	roll 112 mm	roll 58 mm
Leads printing format	1x12, 1x6, 1x4, 1x3, 2x6+1R, 2x6+2R, 4x3+1R, 4x3+2R, 4x3+3R; synchronized/in-time; lead system: standard or Cabrera		2x6 (SD), 4x3, 12x1
PC Connectivity	Cable, WiFi/LAN module (optional)		Cable
Upgrades	Spirometry		N/A

ATTRACTIVE EXTENSIONS



Vacuum electrodes



WiFi/LAN module



Spirometry module



Diagnostic software